FILE NO: 1228.42360X00

CLIENT REF:

INTRAMEDULLARY LOCKED COMPRESSION SCREW FOR STABILIZATION AND

UNION OF COMPLEX ANKLE AND SUBTALAR DEFORMITIES

**TECHNICAL FIELD** 

[0001] The present invention relates generally to the stabilization and reconstruction

of deformities of body joints, and more particularly to double-threaded screws used for

stabilization and reconstruction of deformities of the ankle and subtalar joints, and which

are particularly suitable for stabilization and fusion of the tibiotalar, talocalcaneal,

tibiocalcaneal and/or tibiotalocalcaneal joints.

**BACKGROUND OF INVENTION** 

Severe pain, instability, deformity and/or difficulty walking or bearing weight [0002]

on a joint or joints are usually secondary to painful traumatic arthritis, neuropathic

destruction of joints, post infection arthritis, failed joint replacement, especially of the

ankle and are commonly caused by garden variety primary osteoarthritis or

degenerative arthritis. Such pain, instability, deformity and/or difficulty walking or

bearing weight, typically indicate that fusion in the ankle or subtalar joints is required.

[0003] A wide variety of traditional techniques and devices have been used to effect

fusion in the foot and ankle, such as external fixators, crossed large internal screws

passed from one bone to another, metal plates that pass from one bone to the another

and held in place with screws to each involved bone, or cortical bone grafts, allografts or

FILE NO: 1228.42360X00

CLIENT REF:

autografts that bridge the involved bones and are placed with or without screws fixing

the graft to each involved bone. Generally, to fuse the joints, e.g. the tibiotalar,

talocalcaneal, tibiotalocalcaneal or tibiocalcaneus jioint, as much cartilage as possible is

removed from the adjacent bones of the subject joint. These bones are then coapted or

placed together and then held together tightly.

Many of the traditional techniques for fusion in the ankle and subtalar joints [0004]

require that there be no weight bearing load placed on the involved joint until fusion has

be achieved. This is quite difficult for many patients who have suffered trauma involving

the torso, the spine and the upper extremities, since it makes it very difficult for them to

balance themselves on one lower extremity using a walking aid. This is also quite

difficult, or even impossible, for the frail and the elderly patients.

[0005] Furthermore, surgeons using traditional techniques for reconstruction of

complex painful deformities of the ankle and/or subtalar joints, have experienced a high

percentage of poor results. Often, patients have undergone multiple procedures without

success leaving the soft tissue envelope compromised, and prone to infection and

increased morbidity. In many of the most difficult cases, the residual pain and deformity

has necessitated amputation as the final solution.

[0006] Recently, a technique using intramedullary nails has become popular.

However, conventional nails require the use of elaborate jigs to precisely insert locking

screws through the small holes in the nails, and thereby lock the nail in place. This

FILE NO: 1228.42360X00

CLIENT REF:

makes the insertion of the locking screws very difficult for some surgeons. This is

unfortunate because, although the use of intramedullary nails allows a weight bearing

load to be placed on the involved joint prior to fusion being achieved, the installation

procedure for is much more difficult than the procedure required by traditional

techniques.

[0007] An even more recent alternative to traditional techniques for stabilization

and/or amputation is coaptation across the prepared bones of the joints involved in

fusion, using an intramedullary, double-threaded compression screw. This alternative is

generally used in only very complex and difficult cases of arthropathy of the joints, such

as those mentioned above, which must otherwise be treated by amputation, since

salvage by traditional approaches are usually unsuitable or ineffective.

[8000] The procedure utilized to insert a conventional double-threaded compression

screw into the intramedullary canal is quite simple. An anterior or posterior longitudinal

midline incision is used to access the joint and the cartilaginous surfaces of the joint

involved in fusion. Then a 1.0 inch or so transverse incision is made at the intersection

of a line along the anterior border of the fibula and proceeding across the plantar

surface of the foot, with a line drawn through the center of the plantar aspect of the heel.

Blunt dissection is made down to the inferior surface of the calcaneus. A periosteal

elevator is used to gently push aside the soft tissue from the proposed entry site, in

order to insert a suitable guidewire, which will accommodate the reamer and the

intramedullary double-threaded compression screw.

FILE NO: 1228.42360X00

CLIENT REF:

[0009] The guide wire is inserted from the inferior surface of the calcaneus through

the calcaneus, if a talocalcaneal fusion is planned. Alternatively, the guidewire is

inserted into the intramedullary area of the tibia and across the talus with the foot held in

neutral position of flexion, extension, varus, valgus and rotation, if a tibial-talar fusion is

planned. The intramedullary position of the guide wire is verified by intraoperative x-ray

images or by fluoroscopy, thereby confirming coaptation and alignment.

[0010] The double-threaded compression screw is seated using a cannulated

reamer of an appropriate size. The reamer size is dependent on the size of the

extremity. The reamer is placed over the guide wire to prepare the intramedullary

canals of the bones involved in the fusion. Usually, a cannulated reamer, which is

slightly larger than the reamer used to prepare the tibial canal, is used to prepare the

calcaneal canal entry point. After seating the double-threaded compression screw, a

few threads of the trailing end of the screw will typically remain outside the calcaneus.

The proximal end of the screw must engage the diaphysis of the tibia.

[0011] Rigid fixation is immediate and coaptation of the denuded surfaces is very

precise. Autogenous bone grafting placed at the fusion site is generally recommended

for these complex cases, although not mandatory. Allograft is also acceptable for

placement at the fusion sites.

[0012] However, unlike intramedullar nails used for ankle fusion, this double-

threaded compression screw does not require the use of an elaborate jig to insert one

or more interlocking screws. At present, the intramedullary nails used for ankle fusion,

even with interlocking screws, do not bind well in the calcaneus. Therefore, early

weight bearing with these nails is not wise. Needless to say, insertion of interlocking

screws, even with the elaborate jigs, is difficult for many surgeons. The use of a double

threaded, intramedullary compression screw, because of its firm binding in the

calcaneal cortex compared to a nail without threads, will allow early weight bearing prior

to fusion being achieved. The installation procedure of the double threaded screw for

ankle fusion is much easier than the procedure for conventional intramedullary nails

because the interlocking screws are inserted without the use of a jig and enhance the

fixation already achieved by the threads of the screw.

**OBJECTIVES** 

100131 Accordingly, it is an object of the present invention to provide an improved

lockable double-threaded intramedullary compression screw for stabilizing by fusion, a

body joint such as the ankle and/or subtalar joint.

[0014] It is another object of the present invention to provide a lockable double-

threaded intramedullary compression screw which can be locked in place by one or

more locking screws without the use of an aiming jigs or other alignment device to insert

a locking screw(s).

[0015] It is a further object of the present invention to provide a more easily:

installable implant, which can be used to salvage a lower extremity damaged ankle or

FILE NO: 1228.42360X00

CLIENT REF:

subtalar joints by stabilizing the joints chosen for fusion.

[0016] It is a still further object of the present invention to provide a lockable double-

threaded intramedullary compression screw which can be locked in place to achieve

arthrodesis of the ankle by the insertion of one or more locking screws using palpation

or direct vision.

These and other objects and advantages of the present invention will become [0017]

apparent from the following description of the invention and the appended claims taken

in conjunction with accompanying drawings.

**SUMMARY OF THE INVENTION** 

In accordance with the invention, an intramedullary lockable compression [0018]

screw is provided for stabilizing a joint in a body. The compression screw has an

elongated tubular member and a through-hole configured to accommodate a locking

screw.

[0019] The elongated tubular member extends along a substantially straight first

longitudinal axis between the leading end and the trailing end of the tubular member.

The tubular member includes at least three distinct portions. A threaded leading end

portion has a first diameter and is disposed proximate to the leading end. A threaded

trailing end portion has a second diameter, which is larger than the first diameter, and is

disposed proximate to the trailing end. An unthreaded shaft portion interconnects the

FILE NO: 1228.42360X00

CLIENT REF:

threaded leading end and the threaded trailing end portions.

[0020] - The through-hole is configured to accommodate a locking screw, preferably a

4 mm locking screw, and extends along a straight second longitudinal axis between first

and second openings in the outer periphery of the tubular member. The second

longitudinal axis intersects the first longitudinal axis at an angle of other than 90

degrees. The first opening is in a first area of the outer periphery disposed proximate to

the trailing end of the tubular member. The second opening is in a second area of the

outer periphery disposed distal, as compared to the first opening, to the trailing end of

the tubular member.

[0021] According to an aspect of the invention, the elongated tubular member has a

sleeve portion extending from the trailing end, away from the leading end. The sleeve

portion is aligned with the through-hole and configured to guide the locking screw into

the first opening.

[0022] Typically, the compression screw will additionally include another through-

hole configured to accommodate another locking screw, also preferably of 4 mm. This

other through-hole extends along a straight third longitudinal axis between third and

fourth openings in the outer periphery of the tubular member. The third opening is in a

third area and the fourth opening is in a fourth area of the outer periphery of the tubular

member. Both the third and the fourth areas are proximate to the leading end of the

tubular member. The third longitudinal axis intersects the first longitudinal axis,

FILE NO: 1228.42360X00

CLIENT REF:

preferably at a 90 degree angle.

[0023] Beneficially, a cross section of this other through-hole, taken along the third-

longitudinal axis, has different first and second dimensions. More particularly, the first

dimension, which is in a direction parallel to the first longitudinal axis and is preferably

about 1 inch in length, is larger than the second dimension, which is in a direction

perpendicular to the first longitudinal axis.

[0024] According to still other aspects of the invention, an aperture is provided in the

trailing end of the tubular member. The aperture is configured to accommodate a tool

for applying a torque to the tubular member. Preferably, the aperture includes both a

non-threaded portion configured to accommodate a tool for applying a torque to the

tubular member, and a threaded portion configured to engage a tool for extracting the

tubular member from a body joint.

[0025] Thus, the lockable, double-threaded compression screw described herein

allows early weight bearing even while healing and fusion are taking place. This

reduces the nursing care required in a supervised facility or at home, since it is much

easier for the patient to engage in the activities of daily living, such as going to the toilet

or to a dining area, independently with or without a walking aid.

[0026] An additional benefit of the described locked double threaded compression

screw is that no jigs are required to place the locking screws. This makes the operative

FILE NO: 1228.42360X00

CLIENT REF:

time much shorter and the surgeon's exposure to fluoroscopic radiation during surgery

much less.

Therefore, the subject double-threaded locked intramedullary compression [0027]

screw for fusion is ideal for very complex cases. The procedure for placement of the

locking screw(s) requires much less radiation during surgery than that required for

placement of a locking screw(s) in a conventional intramedullary nail.

[0028] Furthermore, the subject intramedullary double threaded compression screw

can be locked with one or more locking screws, without the jig(s) required to insert the

locking screw(s) in a conventional intramedullary nail. This is because, the proximal

through-hole, near the trailing end of the tubular member, is disposed and orientated so

as to be visible to the surgeon's naked eye, and the distal through-hole, near the

leading end of the tubular member, is in the form of an enlarged slot that is much larger

than the small holes formed in conventional double-threaded compression screws.

Thus, using the compression screw described herein, the locking screws can be applied

proximally and/or distally using direct vision or by palpation, with or without a guidewire.

**BRIEF DESCRIPTION OF THE DRAWINGS** 

Figure 1 is a longitudinal cross-sectional view of a foot and ankle with guide

wire and double-threaded compression screw introduced through the calcaneus and

into the tibia, in accordance with the present invention.

FILE NO: 1228.42360X00

CLIENT REF:

[0030] Figure 2 depicts a lockable double-threaded compression screw according to

the present invention.

[0031] Figure 3 depicts a side view of the double-threaded compression screw

shown in Figure 2.

[0032] Figure 4 depicts a cross-sectional view of the double-threaded compression

screw shown in Figure 3.

[0033] Figure 5 depicts a view of the trailing end of the double-threaded compression

screw shown in Figure 2.

[0034] Figure 6 depicts the trailing end of the double-threaded compression screw

shown in Figure 2, with the proximal locking screw inserted.

PREFERRED EMBODIMENT OF THE INVENTION

[0035] A critical step in the use of the invention is the creation of the intramedullary

canal or cavity across the joints between the bones involved in the fusion. As shown in

Figure 1, the canal extends from the sole of the foot 110, through the calcaneus 112,

and talus 114, and tibia 118. The canal must be reamed to a sufficient size to accept a

guide wire 120 and to accommodate the lockable double-threaded intramedullary

compression screw 130.

FILE NO: 1228.42360X00

CLIENT REF:

[0036] Referring now to Figures 2 and 3, a lockable double-threaded intramedullary

compression screw 200, in accordance with the present invention, is depicted. As

shown, the compression screw 200 has three distinct sections, which will also be

sometimes referred to as portions. One is a narrower, e.g. smaller diameter, threaded

section 215, which is located proximate to the leading end 205 of the compression

screw 200 and has threads 217. The leading end 205 is the end that is first inserted

into the cavity during surgery. Another is a wider, e.g. larger diameter, threaded section

220, which is located proximate to the trailing end 210 of the compression screw 200

and has threads 222. Still another is an unthreaded shaft section 225, which is located

between threaded sections 215 and 220. The unthreaded shaft section 225 is

preferably narrower, e.g. of smaller diameter, than both the threaded sections 215 and

220.

[0037] As shown in Figure 5, the trailing end 210 of the compression screw 200 is

provided with a socket aperture 510 to accommodate a tool (not shown), e.g. an Allren

type wrench, for applying a clockwise torque to the compression screw 200. The

application of the torque will initially cause the threads 217 to engage the bones in, and

in the vicinity, of the joint. As the torque continues to be applied, the threads 222 will

also engage the bones in the vicinity of the joint. This engagement of the threads and

rotation of the compression screw 200 causes the compression screw to be drawn

further into the cavity, passing through the bones in the joint to be fused.

[0038] Ultimately, the compression screw 200, will be drawn fully across the joint

FILE NO: 1228.42360X00

**CLIENT REF:** 

involved in the fusion, with the threads 217 engaged with the bone on the distal side of

the joint and the threads 222 engaged with the bone on the proximal side of the joint, as

shown in Figure 1. The shape and size of the intramedullary canal, after appropriate

reaming, ensures that the lockable double-threaded compression screw 200 fits firmly

into the intramedullary canal when placed therein over a well-placed guide wire, as is

well understood by those skilled in the art.

[0039] The placement of the double-threaded compression screw 200, after proper

reaming, over a guide wire in the intramedullary area of the involved bones, results in

the application of compressive forces across the involved joints, because of the double

threading of the compression screw 200. As described above, the leading threads 217

fix themselves in the intramedullary area of the distal bone involved in the fusion, e.g.

the tibia, fibula and/or talus. The trailing threads 222 pull, coapt and compress the

proximal bone or bones, e.g. the calcaneus and talus, together to form a solid fusion

mass of one bone.

[0040] The addition of locking screws 235 and 245, as shown in Figures 4 and 6,

proximally via through-hole 240 and distally via slotted through-hole 230, prevents the

loss of coaptation and compression. The locking screws 235 and 245 also remove the

influence of any torsion forces. Preferably, the through-holes 230 and 240 are sized to

accommodate a 4 mm or larger locking screw, and accordingly the locking screws 235

and 245 have a threaded diameter at least 4 mm.

FILE NO: 1228.42360X00

CLIENT REF:

The through-hole 240 extends along a straight longitudinal axis between a [0041]

first opening 242, in a first area of an outer periphery of the compression screw 200

proximate to the trailing end portion 220, and a second opening 244, in a second area

of the outer periphery of the compression screw 200 which is further from the trailing

end 210 than the first opening 242. As shown, the second opening 244 is in the

unthreaded shaft portion 225 of the compression screw 200. However, it will be

recognized that in some implementations, the second opening 244 could be located in

the trailing end portion 220 of the compression screw 200. The longitudinal axis of the

through hole 240 intersects the substantially straight longitudinal axis of the

compression screw 200 at an angle of other than 90 degrees. Most typically, this angle

will be 45 degrees or less.

[0042] Extending from the trailing end 210, and away from the leading end of the

compression screw 200, is a sleeve portion 214. The sleeve portion 214 is aligned with

the opening 242 of the through-hole 240 and is designed to guide the locking screw 245

into the first opening 240.

[0043] The through-hole 230 extends along a straight longitudinal axis between a

third opening 232, in a third area of the outer periphery of unthreaded shaft portion 225

of the compression screw 200, and a fourth opening 234; in a fourth area of the outer

periphery of unthreaded shaft portion 225 of the compression screw 200. Both the third

and the fourth areas are proximate to the leading end 205 of the compression screw

200. The longitudinal axis of the through-hole 230 intersects the longitudinal axis of the

FILE NO: 1228.42360X00

CLIENT REF:

compression screw 200. Most typically, this angle will be 90 degrees.

[0044] As shown, a cross section of the through-hole 230, taken along the its

longitudinal axis, has a dimension, preferably of approximately 1 inch, in a direction

parallel to the longitudinal axis of the compression screw 200 which is larger than a

dimension in direction perpendicular to the longitudinal axis of the compression screw.

[0045] Referring to Figures 2-4, leading end portion 215 and the trailing end portion

220 of the lockable double-threaded compression screw 200 are axially spaced apart by

the substantially cylindrical shaft portion 225. As depicted, the leading end portion 215

and the trailing end portion 220 have screw threads 217 and 222 that are similar, but

not necessarily of the same pitch. The leading end 205 of the compression screw 200 is

lodged into the proximal bone involved in the fusion and the threads 222 of the trailing

end portion 220 of the compression screw are, for example, lodged in the inferior area

and cortex of the calcaneus.

[0046] Preferably, the pitch of the threads 217 of the leading end section 215 is

slightly greater than the pitch of the threads 222 of the trailing end section 220 and, as

noted above, the external diameter of the trailing end portion 220 of the compression

screw 200 is greater than the diameter of the leading end portion 215. However, it

should be understood that the relationships between pitch of the threads 217 and 222

and the diameter of the leading and trailing end section 215 and 220 can be varied as

may be desired for the particular implementation.

FILE NO: 1228.42360X00

**CLIENT REF:** 

[0047] As previously discussed, the central area of the trailing end 210 has an

aperture 510, which may be in the form of a slot, disposed in the head 212 of the

compression screw 200. As shown, the head 212 forms part of the trailing end portion

220 of the compression screw 200. The aperture 510 is configured to accommodate an

insertion tool (not shown), such as a wrench or screwdriver, for applying a clockwise

torque to the compression screw 200. Application of the clockwise torques drives the

compression screw 200 into the bone or bones involved in the fusion.

[0048] The aperture 510 may also accommodate an extraction tool (not shown),

such as a wrench or screwdriver, for applying a counter-clockwise torque to the

compression screw 200. The extraction tool and insertion tool may be the same tool or

different tools. The application of the counter-clockwise torque drives the compression

screw 200 from its position in the bone or bones involved in the fusion.

[0049] The compression screw 200 illustrated in Figures 1-6 is principally to be

lodged in spongy inner tissue of a bone. The helical threads 217 and 222 create

compression between the bones involved in the fusion.

[0050] Optionally, internal threads 250 are provided beyond the aperture 510 in the

head 212 of the compression screw 200. The internal threads 250 are formed within

the trailing end portion 220 of the compression screw 200. These threads are adapted

to accommodate a threaded extraction tool (not shown). The threads of the threaded

extraction tool will engage the internal threads 250 as the threaded extraction tool is

FILE NO: 1228.42360X00

CLIENT REF:

turned counter-clockwise, thereby locking the threaded extraction tool to the

compression screw 200. After locking, the threaded extraction tool can be used to

apply a further a counter-clockwise torque to the compression screw 200.. The

application of the counter-clockwise torque, after locking, drives the compression screw

200 from its position in the bone or bones involved in the fusion. Furthermore, because

the threads of the threaded extraction tool are in locked engagement with the internal

threads 250, the threaded extraction tool can also be used to apply a tension force to

the compression screw 200. The application of the tension force, after locking, may be

sufficient to pull the compression screw 200 from its position in the bone or bones

involved in the fusion.

[0051] The use of the compression screw 200, and in particular the extraction of the

compression screw 200 from the bones involved in the fusion, is greatly facilitated by

the accommodation of an extraction tool that can be locked into the trailing end portion

220 of the compression screw 200 by counter-clockwise engagement with the internal

threads 250, via the trailing end 210 of the compression screw.

[0052] In using the double-threaded compression screw 200, it is first necessary to

select a compression screw suited to the task at hand, and to then follow the insertion

procedure described above. The selection of the screw can be done by measuring

preoperative x-ray images of the bones involved.

[0053] The foregoing description of the present invention is illustrative and

DOCKET NO: 3023-004 FILE NO: 1228.42360X00

CLIENT REF:

explanatory. Various changes in size, shape and materials, as well as in the details of

the preferred embodiment of the illustrated construction thereof may be made without

departing from the spirit of the invention.

[0054] Referring again to Figure 1, an anterior longitudinal midline approach to the

ankle was used to approach the joint and to débride the cartilage from the distal tibia

and the proximal talus. Next, a 1 inch transverse plantar incision was made at the

intersection of a line drawn along the anterior border of the fibula proceeding along the

plantar surface of the foot with a line drawn to the plantar center of the heel. Blunt

dissection was made down to the inferior surface of the calcaneus. A periosteal elevator

was used to gently push the soft tissue from the proposed entry site for an appropriate

3/32 guide wire to be centered through the intersection process of the heel and through

the calcaneus, through the talus and into the intramedullary area of the distal tibia.

While the guide wire was inserted, the foot was held in neutral position of flexion,

extension, varus, valgus and rotation. The intramedullary position of the guide wire was

verified by intraoperative xray images or fluoroscopy, and coaptation and alignment

confirmed. A 7 mm reamer was used to prepare the path for the double-threaded

compression screw being inserted over the guide wire. A 10 or 11 mm reamer was used

to prepare the calcaneus. A compression screw of sufficient length, was selected,

based on preoperative x-ray images of the involved bone. In this case the selected

compression screw had a length of 6 inches. The double threaded screw must be long

enough to engage the diaphysis of the tibia. The selected compression screw was

driven over the guide wire through the calcaneus and into the intramedullary space of

FILE NO: 1228.42360X00

CLIENT REF:

the distal tibia. Coaptation of the prepared surfaces of the distal tibia and proximal talus

was rigid and immediate. Allograft was placed at the fusion site.

[0055] Using direct vision, a locking screw was passed over a guide wire into the

proximal through hole, in the trailing end portion of the double-threaded compression

screw. By palpation of the distal slot proximate to the leading end portion of the

compression screw, a guide wire was inserted. The position of this guide wire passing

through the distal slot in the leading end portion of the double-threaded screw was

confirmed anteriorly, posteriorly and laterally by fluoroscopy. The length of the guide

wire was confirmed by measurement of the exterior length of the guide wire using a

similar length guide wire. A locking screw of appropriate length was then selected and

inserted over the guide wire. A torque was applied to the locking screw to thereby

securely distally lock the compression screw in position. The wounds were irrigated

thoroughly, more allograft bone was then placed around the fusion site, and the wounds

were closed in the usual fashion. A sterile bulky dressing was applied and a well-

padded short leg cast was then applied. Partial weight bearing on the lower extremity,

using a walking aid, was allowed immediately.